The importance of sonography

Treatment of diseases of the musculoskeletal system with **ESWT**

Matthias Beck

Targeted and precise application of focused extracorporeal shock waves in the treatment of some ESWT indications (calcified tendinitis, heel spur, patellar tendinitis) is difficult or impossible without the use of sonography. While B-mode sonography used to focus on the detection and localisation of tissue alterations, the implementation of power Doppler sonography now enables identification of changes in microcirculation and perfusion associated with chronic inflammations in soft tissue (tendons, tendon sheaths, joint capsules, ligaments) (Wakefield RJ 2003, Belli P 2000, Wakefield R 2000, Klauser AE 2002). Contrary to colour flow Doppler (CFD) imaging, power Doppler sonography (PDS) does not represent the calculated mean value of the Doppler frequency shift (determined by the direction and velocity of the blood cells), but rather displays the amplitude of the Doppler signal as a colour spectrum, providing information on the blood volume present in the tissue at the time of imaging. Consequently, CFD is indicated for the

examination of larger vessels (e.g. carotid artery), while PDS can better display slower flow velocities in smaller vessels (e.g. in joint capsules and tendon sheaths). Especially in the examination of enthesopathies, D'Agostino et al. (2003) were able to identify a positive correlation between flow changes (pathological vascularisation) in the junction region of the tendon with the cortical bone tissue and the degree of tendon disease. Based on their observations, they developed a 5-stage classification system of enthesopathies.

Fig.: Sonographic examination of Achilles tendon with the new DUOLITH® SD1 »ultra« from STORZ MEDICAL

This represents a major benefit for validating the diagnosis or differential diagnosis of diseases of the musculoskeletal system and for confirming the regression of inflammation-induced alterations during ESWT as an indication of tissue improvement achieved with this therapy (Walther M 2001, Hau M 1999, M. I. Richards PJ 2009, D. A. Richards PJ 2001, Hoksrud A 2008).

Studies so far conducted on this issue have revealed a high predictive value for Achilles tendon inflammation (M. I. Richards PJ 2009, D. A. Richards PJ 2001, Reiter M 2004) and patellar tendinitis (Hau M 1999, Hoksrud A 2008). Other indications such as radial tunnel syndrome (Bodner G 2002), inflammation of sacroiliac joints (Klauser A 2005) and radial epicondylopathy (Clarke AW 2010) have been investigated in isolated studies, which, however, have a low evidence class so far. The relatively new method of contrast medium enhancement in the tissue to be examined will certainly further improve the validity of sonography in the years to come (Hau M 1999, D. A. Richards PJ 2001, Hoksrud A 2008) and thus optimise the existing link between sonography and ESWT.

Thanks to the new generation of STORZ MEDICAL shock wave therapy systems, the technical requirements for this link are already in place.

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Source: STORZ MEDICAL AG

To view the full version of this article and the detailed reference list (in German), visit <http://bit.ly/H8u5VV> or www.ortho-online.de.

Fig.: PDS of an inflammatory paratenon change of the Achilles tendon in a marathon runner

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GmbH & Co. KG
► Quelle: STORZ MEDICAL AG