Getting the application pressure right

Radial shock wave therapy in orthopaedics

Extracorporeal shock wave therapy (ESWT) was first used in pain therapy in the mid-1990s. However, owing to their shape and cost, the urological lithotripters used for this purpose did not manage to work their way into medical practice. In 1999, low-energy radial shock waves (R-SW) started their success story in the field of orthopaedics with compact, conveniently-priced and practice-oriented shock wave systems.

Radial shock waves (R-SW) are applied to soft tissue in a gentle manner and without any lasting side effects. The success of the treatment depends on the application pressure (1.4 – 5 bar), the coupling transmitter employed and on the shock wave frequency (1 – 21 Hz). Experience gathered so far has shown that an application pressure below 1.4 bar will not provide sufficient effectiveness of the therapy. On the other hand, application pressures in excess of 4 bar have shown not to improve treatment results any further. Also, shock frequencies below 8 Hz are not advised.

Fig. 1: R-SW handpiece with innovative handpiece display for pressure and frequency selection

The different materials employed to couple the shock wave energy to the patient's body have been found to produce different effects. If a ceramic coupling medium is used (C15 CERAmA-x[®] transmitter), for example, patients perceive less pain during treatment. Titanium (D20-T transmitter) allows the energy level and penetration depth to be increased, while transmitters consisting of a focused lens filled with special fluid (F15 focus lens transmitter) are ideal for the gentle treatment of near-surface body regions (e.g. face).

Fig. 2: R-SW transmitters (from left to right): R15 ESWT transmitter (Ø 15 mm), C15 CERAmA-x[®] transmitter (Ø 15 mm), D20-T transmitter (Ø 20 mm), F15 focus lens transmitter (Ø 15 mm)

Typical R-SW indications include tennis elbows and heel spurs, for example. In addition to their instant analgesic effect during the therapy session, radial shock waves have shown to provide good or even excellent long-term results.

Fig. 3: Treatment of shoulder with R-SW

In 2005, the valid and safe procedure of R-SW therapy was extended to the treatment of myofascial pain syndromes. This was done by adapting the therapy systems to more demanding requirements and by paving the way for new therapy approaches with new technology solutions. As a result, the mere treatment of individual painful spots in tendon insertions can now be complemented by large-area vibratory muscle relaxation therapy. The evolution of R-SW technology has been strongly shaped by specific user experience, which has been consistently translated into technological solutions. This experience has been corroborated by scientific studies. Each patient perceives pain differently and responds differently to shock wave treatment. As a result, SW therapy today is no longer used as an isolated treatment modality but forms part of a complex, individual and holistic therapy concept.

Correct diagnosis and differential diagnosis are essential requirements for successful treatment. Ultrasound and, where necessary, X-ray examinations should be performed.

Fig. 4: DUOLITH® SD1 »ultra« with integrated ultrasound diagnosis, ultrasound image (shoulder front view)

Moreover, the documentation of the therapy process is gaining increasing importance for patients and health insurance providers. To accommodate this requirement, state-of-the-art therapy systems are equipped with on-board patient documentation solutions. Hardly any other therapy approach has developed so dynamically and at such a pace in only twelve years. What needs to be done now is to implement and further develop the experience gathered from observational and scientific studies even more effectively and to avoid placebo effects or poor treatment results which may result from low penetration depths and energy bundling or from excessive application pressures (over 4 bar). In the end, what really matters is the therapy outcome for the patient, who, after all, completes the circle between users and industry.

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(Chiemgau Medical Services – Orthopaedics & Surgery)

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► **Quelle:** STORZ MEDICAL AG