

IDS CRYO-ZET

AIR COOLING SYSTEM



ids

Innovation for better life

The most powerful
Air Cooling System

NEW REVOLUTION!

01

Convenient and economical for users

- It is very economical because treatment times can be considerably reduced.
- It is very convenient because using is very easy.
- It is very useful because consumable parts are not needed.

02

Very effective and safe for patients

- It can reduce skin damage because it can be used with lasers.
- It can reduce pain and recovery time.
- It can protect skin because cooling can be made evenly.

03

The suitable lasers

- Fractional CO2 Laser, IPL, RF Microneedle type, Er YAG, Diode laser, Alexandrite, Nd:YAG(Q-Switched, Long Pulsed YAG)

04

Special Features

- ✓ Air flow controlling system.
- ✓ Cooling air temperature (Max. -35 degree).
- ✓ Quick cooling with high quality compressor.
- ✓ Auto adjustable controlling system.
- ✓ Auto detecting system to keep best condition.

Cryo-Zet



Technical Specifications

A.C Input	220~240V, 50/60Hz	Air Flow Output	1000L/Min
Air Level	10 Steps	Temperature	0~-35°C
Size	400(W)x610(L)x963(H)	Weight	50kg