


GaitKeeper GK-S20 / S22



	GK-S20	GK-S22
Walking Surface Width	20	22
Walking Surface Length	56	56
Treadmill Width	27.5	29.5
Treadmill Length	73	73
Maximum User Weight	400 lbs	400 lbs*
Motor	AC 2HP	AC 2HP
Incline Range	0-15%	0-15%
Forward Speed Range (mph)	0.1 - 6.0	0.1 - 10.0
Reverse Speed Range (mph)	0.1 - 3.0	0.1 - 3.0
Footprint	73 x 32 x 53	73 x 33 x 52
Step-Up Height	7.5	7.5
Handrail Length	37	37
Handrail Height (MAX)	40	40
AdjustaBars Range	24-40	24-40
GaitSens	N/A	Option
Shipping Dimension	76 x 35 x 19	77 x 36 x 20
Shipping Weight	312 lbs	324 lbs

*can be upgraded to 500 lbs

All dimensions subject to change.
Dimensions listed in inches

Features:

- True Zero Start
- Removable Console
- AdjustaBars™
- Reverse Feature

Power Requirements	110 V - 120VAC, 60Hz, 15 amp
Regulatory Approvals	TUVc, TUVus, Standard UL, CAN/CSA,



GK-S22 w/mounted
GaitSens Tablet

