


GaitKeeper GK-S20



	GK-S20
Walking Surface Width	20
Walking Surface Length	56
Treadmill Width	27.5
Treadmill Length	73
Maximum User Weight	400 lbs
Motor	AC 2HP
Incline Range	0-15%
Forward Speed Range (mph)	0.1 - 6.0
Reverse Speed Range (mph)	0.1 - 3.0
Footprint	73 x 32 x 53
Step-Up Height	7.5
Handrail Length	37
Handrail Height (MAX)	40
AdjustaBars Range	24-40
GaitSens	N/A
Shipping Dimension	76 x 35 x 19
Shipping Weight	312 lb

All dimensions subject to change.
Dimensions listed in inches

Features:

- True Zero Start
- Removable Console
- AdjustaBars™
- Reverse Feature

Power Requirements	110 V - 120VAC, 60Hz, 15 amp
Regulatory Approvals	TUVc, TUVus, Standard UL, CAN/CSA,

