


GaitKeeper GK-S22



GK-S22

Walking Surface Width	22
Walking Surface Length	56
Treadmill Width	29.5
Treadmill Length	73
Maximum User Weight	400 lbs*
Motor	AC 2HP
Incline Range	0-15%
Forward Speed Range (mph)	0.1 - 10.0
Reverse Speed Range (mph)	0.1 - 3.0
Footprint	73 x 33 x 52
Step-Up Height	7.5
Handrail Length	37
Handrail Height (MAX)	40
AdjustaBars Range	24-40
GaitSens	Option
Shipping Dimension	77 x 36 x 20
Shipping Weight	324 lb

*can be upgraded to 500 lbs

All dimensions subject to change.
Dimensions listed in inches

Features:

- True Zero Start
- Removable Console
- Adjustabars™
- Reverse Feature

Power Requirements	110 V - 120VAC, 60Hz, 15 amp
Regulatory Approvals	TUVc, TUVus, Standard UL, CAN/CSA,



GK-S22 w/mounted
GaitSens Tablet

