

GaitKeeper GKmini

Features:

- Portable Design for use with LiteGait®
- Hand Held Remote with LED Display
- High torque drive motor with 0.1 mph starting speed with 0.1mph increments



GaitKeeper Mini 21	
Walking surface width	21
Walking surface length	44
Treadmill width	24
Treadmill length	47
Footprint	47 L x 24 W
Maximum user weight	350
Handheld remote	Yes
Handlebars	Removable
Handlebar height	37
Incline range	N/A
Drive motor	Brushless
Speed range	0.1 – 3mph
Step-up height	6

Accessories

GaitSens 2.0

Optional

Power Requirements	110-240 VAC 6A, 50/60 Hz
Regulatory Approvals	CE-60601
Shipping Weight	110

All dimensions subject to change.
Dimensions listed in ./lbs./mph.

GaitKeeper™ is owned by

mobility research



Mobility Research Inc

LiteGait.com · 3156 w. lewis avenue · phoenix 85009, az
ph. 800.332.9255 or 480.829.1727 · fax 480.829.0737 · sales@LiteGait.com